

RESPIRATORY COMPROMISE AND CAPNOGRAPHY

WWW.PROMISETOAMANDA.ORG/AMANDAS-STORY/

Medtronic

RESPIRATORY COMPROMISE¹

“Respiratory Compromise (RC) is a potentially life threatening state of unstable respiratory health: It is a gradual, subtle imbalance in patient response that encompasses respiratory failure and arrest, with symptoms that manifest differently in each patient.”

-Respiratory Compromise Institute
www.respiratorycompromise.org

Medtronic

RESPIRATORY COMPROMISE CATCHING THE CASCADE

01 **At Risk**
Patient Desaturates

02 **Respiratory Insufficiency**
2LPM NC, continues to desaturate
4LPM NC, continues to desaturate

03 **Respiratory Failure**
6LPM NC, continues to desaturate
NRB applied

04 **Respiratory Arrest**
Patient resting comfortably in nurse notes

05 **Airway/Vent Support Needed. Could be found "Dead in Bed"**

Medtronic

RESPIRATORY COMPROMISE BY THE NUMBERS

SECOND leading avoidable patient safety issue

THIRD most rapidly increasing hospital inpatient cost in the US

THIRD most common cause of avoidable deaths in the U.S. (estimated 300,000-400,000 per year)

FIFTH top condition leading to increasing hospital costs¹

Risk of death is **29%** higher among patients with respiratory compromise

Common.
Costly.
Deadly.
Preventable.

Medtronic

CAPNOGRAPHY

Medtronic

OVERVIEW

What is Capnography

Capnography is the non-invasive, continuous measurement of CO₂ concentration at the airway

Provides:

- Respiratory rate detected from the actual airflow
- Numeric "end tidal carbon dioxide" (etCO₂)
 - Normal text book range 35-45 mmHg
- Waveform Representing CO₂ concentration in every breath

Medtronic

Obtaining an Accurate Respiratory Rate

Manual Counting

Measures:

- Chest or air movement
- Based on observation or auscultation that may be restricted by patient movement, draping or technique

Impedance (ECG Leads)

Measures:

- Attempt to breathe
- Chest movement
- Based on measuring respiratory effort or any other sufficient movement of the chest

etCO₂

Measures:

- Actual exhaled breath at airway
- Hypoventilation and No Breath detected immediately!
- Most accurate RR, even when you are not in the room!

Medtronic

KNOWLEDGE TEST 1

All of the following are patient associated risk factors for RC following opioid administration except:

- Patient with BMI of > 30kg/m²
- Central sleep apnea patient
- Patient with kyphoscoliosis
- Intraoperative administration of ketamine ←
- Use of PCA on a patient
- Patient with pulmonary sarcoidosis
- CHF patient

Medtronic

RESPIRATORY CYCLE

The process of getting O₂ into the body

Oxygenation

Z- separate physiologic processes

SpO₂

etCO₂

Ventilation

The process of eliminating CO₂ from the body

Medtronic

DIFFERENTIATION OF SpO₂ & EtCO₂

Pulse Oximetry

Measures SpO₂ is saturation level in the blood

- Values lag with hypoventilation & apnea
- Supplemental O₂ increases this lag time

Capnography

Measures etCO₂ ventilation

- ventilation & hypoventilation & apnea detected immediately
- Respiratory Depression is 28 times more likely to be detected with Capnography than with Pulse Oximetry*

*Measuring to Improve Ventilation: Safety, Quality, Satisfaction and Alignment. American Society of Anesthesiologists (ASA). 2014. ASA. www.asahq.org/ASA/education/ASA-CO2-Check. ©ASA. All rights reserved and/or assigned, 2014.

Medtronic

CAPNOGRAPHY ASSESSES PHYSIOLOGICAL FUNCTIONS

METABOLISM

Determines CO₂ production

PERFUSION

How efficiently CO₂ is returned to the lungs

VENTILATION

Adequacy of lung function (COPD, stiff lungs, fluid or sticky lungs)

ALL THREE ARE IMPORTANT!

Medtronic

FACTORS THAT INCREASE ETCO2


Metabolic	Perfusion	Ventilation	Technical Errors
Fever Malignant hyperpyrexia Sodium bicarbonate Tourniquet release Venous CO ₂ embolism	Increased cardiac output Increased blood pressure ROSC Effective chest compressions	Bronchial intubation Partial airway obstruction Rebreathing (expiratory time too short) Hypoventilation w/sedation RR or Vt is too low on ventilator	Exhausted CO ₂ absorber Inadequate fresh gas flows Leaks in breathing system Faulty ventilator Faulty valves

Medtronic

FACTORS THAT DECREASE ETCO2

Metabolic	Perfusion	Ventilation	Technical Errors
Hypothermia NPO	Reduced cardiac output Hypotension Hypovolemia Pulmonary embolism Cardiac arrest	Shallow breathing Hyperventilation Apnea Total/Partial airway obstruction Tracheal extubation V/Q mismatch or Shunting/ARDS RR or Vt too high on MV	Circuit disconnection Sampling tube leak or kinked ET tube Malfunction of ventilator

Medtronic



VENTILATION & PERFUSION MATCHING

Normal Shunt Perfusion

- Physiologic Dead Space
- Normal Physiologic Shunt

(a) EtPCO₂ = 2.5 mmHg

Deadspace Ventilation

- Ventilation - No Perfusion
- EtCO₂ < PaCO₂

- Pulmonary embolism
- Cardiac arrest
- Hypovolemia
- Hypoxic Pulmonary Vasoconstriction

Abnormal Shunt Perfusion

- Perfusion - No Ventilation
- EtCO₂ < PaCO₂

- Pulmonary Edema
- Pneumonia
- ALI/ARDS
- Bronchial Intubation
- Mucus Plugging
- Bronchospasm
- Severe Atelectasis

Medtronic

Name: Leah Coufal
Age: 11
Diagnosis: Pectus carinatum repair

"Would real-time monitoring have saved Leah? That is one of the many questions that I have asked myself every day since I found my daughter, Leah, dead in her hospital bed. The answer is yes, it would have." - Lenora Alexander

http://www.medscape.com/viewarticle/906000_1
<http://www.scribd.com/doc/100000000/Leah-Coufal>

Medtronic

MANY CONDITIONS MAY CAUSE RESPIRATORY FAILURE...


Questions/Complexity arise

- What do I monitor COPD vs pneumonia vs postop vs HF vs...?
- How does monitoring guide intervention?

Medtronic

HYPOXEMIC RESPIRATORY FAILURE (TYPE I) OXYGENATION FAILURE

- Primary Issue** – Decreasing ability to oxygenate blood due to congestion (e.g., fluid, pus, etc.) in lung tissue (alveoli)
- Examples** – Pneumonia, ARDS, sepsis, HF, etc
- Monitor Trend Results**
 - S_pO_2 is maintained above 90% through a significant drop in arterial oxygen (P_aO_2) due to:
 - To compensate, patient breathes faster and/or deeper increasing minute volume (MV) and clinicians give supplemental O_2
 - Decreasing blood CO_2 ($PaCO_2$ or $etCO_2$)
- Interventions** – Aimed at clearing congestion and improving oxygenation (FRC)



Graphs adapted from: Lynn LA, Curry J.P. Patterns of unexpected in-hospital deaths: a root cause analysis. Safety in Surgery 2011, 5:3 (epub) <http://www.pssjournal.com/content/5/1/3>

Curry and Jungquist Patient Safety in Surgery 2014, 8:29 <http://www.pssjournal.com/content/8/1/29>

Medtronic

EXAMPLE OF HYPOXEMIC RESPIRATORY FAILURE (TYPE I) OXYGENATION FAILURE


“A healthy male who had just undergone elective surgery develops shortness of breath that’s noticed by his family who express concern to the nurse. The nurse, citing a normal oxygen saturation reading on his oximeter, reassures the family that the monitor indicates he’s okay. Eventually his respiratory rate does rise to a critical value, but by this time it’s too late to effectively respond to his rapidly deteriorating clinical condition and the patient, with sepsis, dies.”

Curry and Jungquist Patient Safety in Surgery 2014, 8:29 <http://www.pssjournal.com/content/8/1/29>

Medtronic

HYPERCAPNIC RESPIRATORY FAILURE (TYPE II) VENTILATION FAILURE

- Primary Issue**: Inadequate ventilation to clear carbon dioxide – Hypoventilation
- Decreased RR and/or tidal volume = decreased minute volume (MV). RR can be normal and MV is decreasing as a result of shallow breathing
- Examples** – Opioid use, sedation, brain injury, obesity hypoventilation syndrome, etc.
- Monitor Trend Results**
 - S_pO_2 is maintained above 90% through a significant rise in CO_2
 - Normal, decreasing, or increasing RR
 - Increasing CO_2 (ABG or $etCO_2$)
- Interventions** aimed at improving ventilation
 - Airway maneuvers
 - Reducing respiratory depression
 - Ventilatory support



Graphs adapted from: Lynn LA, Curry J.P. Patterns of unexpected in-hospital deaths: a root cause analysis. Safety in Surgery 2011, 5:3 (epub) <http://www.pssjournal.com/content/5/1/3>

Curry and Jungquist Patient Safety in Surgery 2014, 8:29 <http://www.pssjournal.com/content/8/1/29>

Medtronic

EXAMPLE OF HYPERCAPNIC FAILURE (TYPE II) VENTILATORY FAILURE


“A healthy female who is receiving routine post-op nasal oxygen has been up all night complaining of severe postop pain, but is now finally asleep after yet another dose of IV opioid. The nurse, noticing on rounds the patient’s oxygen saturation is ‘perfect’ on the monitor, decides not to awaken her. She is found dead in bed 4 hours later.”

Curry and Jungquist Patient Safety in Surgery 2014, 8:29 <http://www.pssjournal.com/content/8/1/29>

Medtronic

AROUSAL FAILURE (TYPE III) SLEEP DISORDERED BREATHING

- Primary Issue** – Recurring apneas and/or shallow breathing during sleep
- Repetitive pattern of hypoventilation followed by hyperventilation
- Normal protective mechanism is to ‘arouse’ when O_2/CO_2 levels are abnormal
- Recovery response may be blunted by opioids, sedatives, alcohol, etc.
- Examples** – Obstructive sleep apnea, central sleep apnea, Cheyne-Stokes Resp
- Monitor Results**
 - Sawtooth pattern of drops and increases in RR, MV, SpO_2 with reciprocal change in CO_2
- Interventions**, Unidx SDB management
 - Screening tools
 - Positive Airway Pressure therapies
 - Judicious use of arousal/RD suppressants



Graphs adapted from: Lynn LA, Curry J.P. Patterns of unexpected in-hospital deaths: a root cause analysis. Safety in Surgery 2011, 5:3 (epub) <http://www.pssjournal.com/content/5/1/3>

Curry and Jungquist Patient Safety in Surgery 2014, 8:29 <http://www.pssjournal.com/content/8/1/29>

Medtronic

EXAMPLE OF AROUSAL FAILURE (TYPE III) AROUSAL FAILURE

“An otherwise healthy male with unrecognized sleep apnea receives a post-operative opioid. His alarm sounds repeatedly but lasts only for about 30 seconds before it stops, only to repeat again and again. When the nurse awakens the patient he feels fine and is completely alert, asking for more pain medication, which the nurse gives in a normal dose. The nurse, suffering from alarm fatigue, stops responding to the same alarming. Later that night the patient is found dead in bed.”

Curry and Jungquist Patient Safety in Surgery 2014, 8:29 <http://www.pssjournal.com/content/8/1/29>

Medtronic

KNOWLEDGE TEST 2

Hypercapnia can occur with normal SpO₂?

- a) YES
- b) NO




Medtronic

Capnography Waveform 101

Medtronic

NORMAL CAPNOGRAPHIC WAVEFORM



- I: Baseline = no CO₂ in breath
- II: Rapid rise in CO₂ (clearing dead space)
- III: Alveolar plateau (gas exchange)
- IV: End expiration (etCO₂ measured)
- V: Inhalation

Medtronic


NORMAL WAVEFORM



- SpO₂ – Normal
- EtCO₂ - Normal Range 35-45
- Waveform – Normal
- RR – Normal

Medtronic


HYPOVENTILATION WAVEFORM



- SpO₂ – Normal
- EtCO₂ – Increased
- Waveform – increased in amplitude and width
- RR – very decreased

Medtronic

HYPERVENTILATION WAVEFORM



- SpO₂ – Normal
- EtCO₂ – decreasing
- Waveform – decreasing in amplitude and width
- RR – increased

Medtronic

APNEA/LOSS OF WAVEFORM

SpO₂ – Normal or sudden decrease noted
 EtCO₂ – Zero
 Waveform – Absent
 RR – Zero
 Other – no chest wall movement or breath sounds noted

Medtronic

AIRWAY OBSTRUCTION

Absent alveolar plateau indicates incomplete alveolar emptying or loss of airway integrity.
 -Tongue or position of head
 -Secretions
 -Talking or snoring

Medtronic

SHALLOW BREATHING

Hypoventilation with shallow respirations

Medtronic

HYPOPNEIC HYPOVENTILATION (SHALLOW BREATHING) A MOST MISUNDERSTOOD WAVEFORM

- During very shallow breathing, the waveforms don't reach the alveolar plateau (see waveforms A)
- These short, incomplete waveforms can continue for extended periods (several minutes) causing a build up of CO₂ in the blood while etCO₂ is reading low
- Eventually, the patient does take a normal sized breath with an alveolar plateau (B in graphic), the etCO₂ reading will be high and waveform will be tall. The extent of hypercapnia will depend on the length and severity of the shallow breathing.
- FYI: The period of shallow breathing would have to be much longer than just a few breaths to show a low etCO₂ reading and result in hypercapnia.

Medtronic

KNOWLEDGE TEST 3

Cases of respiratory depression are _____ times as likely to be detected, if they are monitored with capnography, as those who were not monitored.

a) 5
 b) 17
 c) 28
 d) 36

←

Medtronic

ABNORMAL WAVEFORM - WHAT TO DO NON-INTUBATED PATIENT

- 01 Assess Patient
- 02 Check FilterLine® position – reposition as necessary
- 03 Check patient's head and neck position – reposition as necessary
- 04 Periodically instruct patient to take a deep breath as necessary
- 05 If patient is "not breathing" Follow protocol (ABC's) CO₂

Medtronic

What Do The Alarms Mean?

Does the Patient have:	Clinical Assessment Considerations:
<ul style="list-style-type: none"> Low RR alarms Decreasing RR 	<ul style="list-style-type: none"> Hypoventilation Shallow breathing Periods of "No Breath" or Apnea Check if etCO₂ level is increasing
<ul style="list-style-type: none"> High etCO₂ alarm Increasing etCO₂ 	<ul style="list-style-type: none"> CO₂ accumulating Low RR Recent ingestion of carbonated beverages/antacids
<ul style="list-style-type: none"> Low etCO₂ alarm Decreasing etCO₂ 	<ul style="list-style-type: none"> RR increasing/hyperventilating Shallow breathing
No Breath or Apnea alarm	<ul style="list-style-type: none"> Apnea Complete airway obstruction FilterLine needs repositioning Extreme shallow breathing

Medtronic

POSSIBLE CLINICAL INTERVENTION

- Ensure an open airway Allows ventilation, thus the removal of CO₂
- Stimulate your patient to take deep breaths Removes CO₂
- Reposition FilterLine® if necessary Provides accurate sample
- Follow your policy Know when to assist ventilation or withhold PCA medications


Remember: CO2 removal must keep-up with production!

Medtronic

SOME KEYS TO SUCCESS


Educate the patient and family before you begin-get buy-in.

- What: "This FilterLine will help monitor your breathing".
- Why: "Your medication may slow your breathing. The cannula allows us to monitor your breathing even when I am not here. If you hear a beep, it is a reminder to take a breath".



It is all about patient safety!

Medtronic



2015 AHA GUIDELINES FOR CARDIOPULMONARY RESUSCITATION AND EMERGENCY CARDIOVASCULAR CARE SCIENCE

"Continuous waveform capnography is recommended in addition to clinical assessment as the most reliable method of confirming and monitoring correct placement of an endotracheal tube."

"...It is reasonable to consider using quantitative waveform capnography in intubated patients to monitor CPR quality, optimize chest compressions, and detect ROSC during chest compressions or when rhythm check reveals an organized rhythm. If PetCO₂ is <10 mm Hg, it is reasonable to consider trying to improve CPR quality by optimizing chest compression parameters. If PetCO₂ abruptly increases to a normal value (35 to 40 mm Hg), it is reasonable to consider that this is an indicator of ROSC."

Medtronic

THE WALL STREET JOURNAL

"If we didn't have the CO₂ readings we were getting, we would have terminated efforts much sooner!" - Bruce Goodman, Mayo Flight Paramedic

Name: Howard Snitzer
Age: 54
Diagnosis: Cardiac Arrest

96 Minutes - On the 12th defibrillation shock after start, a rise in capnography indicates ROSC.



1. Winslow, R. "96 minutes without a heartbeat". *Wall Street Journal*, May 11, 2011.

Medtronic

THE JOINT COMMISSION *SENTINEL EVENT ALERT* ISSUE 49 AUGUST 2012

Regarding the safe use of opioids in hospitals

"In addition to monitoring respiration and sedation, pulse oximetry can be used to monitor oxygenation, and capnography can be used to monitor ventilation. Staff should be educated not to rely on pulse oximetry alone because pulse oximetry can suggest adequate oxygen saturation in patients who are actively experiencing respiratory depression, especially when supplemental oxygen is being used – thus the value of using capnography to monitor ventilation. When pulse oximetry or capnography is used, it should be used continuously rather than intermittently."

Medtronic

USING CAPNOGRAPHY WITH PATIENT CONTROLLED ANALGESIA

Respiratory depression is more likely to occur when stimulus to the patient (i.e., clinicians, family etc.) stops

Accurately monitors effective ventilation and respiratory rate, even when clinicians are away from the room

Earliest indicator of:

- Apnea
- Obstruction
- Hypoventilation

Watch for: No breath (apnea), shallow breathing, ↓ RR & ↑ etCO₂ trend

Titrate dosage (duration between doses or actual dose amount) according to trends

Medtronic

UNIQUE APPLICATIONS OF ETCO₂ IN CRITICAL CARE

- Gastric tube placement
- Reduce ABG's and trend a-A gap for improvement
- Passive leg raise test (sepsis)
- DKA

Medtronic

ETCO₂ PREDICTS FLUID RESPONSIVENESS IN PASSIVE LEG RAISING

- 65 ventilated patients needing volume expansion
- Compared changes in EtCO₂ with arterial pressure to reflect changes in CO
- EtCO₂ increase ≥ 5% predicted fluid responsiveness (p=0.0001)
 - Increase in the CI ≥ 15%
 - Sensitivity 71% (95% CI 48 – 89%) and specificity of 100% (CI 82 – 100%)
- The changes in EtCO₂ induced by a PLR test predicted fluid responsiveness with reliability, while the *changes in arterial pulse pressure did not.*

Medtronic

PASSIVE LEG RAISING

- Patient – HOB 45 degrees
- Obtain Capnography reading
- Initial EtCO₂: 34 mmHg
- Lift legs for 1 – 2 minute
- 40 mmHg after 90 sec
- What does the patient need?
- **Fluids!!!!**

Medtronic

CONCLUSION

Respiratory Compromise is **COMMON, COSTLY, DEADLY & PREVENTABLE**

Capnography provides an earlier detection of Respiratory Compromise

You have the power to save lives!

Medtronic

QUESTIONS AND COMMENTS
THANK YOU!!!