

Polymerase Chain Reaction (PCR): Amplification of insert

- 1) Resuspend primers in buffer EB (elution buffer from Qiagen kit) to a final concentration of 60 μM .
For 25 nmol of primers, 400 μl EB.
- 2) Taq Readymix + MgCl_2 from Sigma contains dNTPs, reaction buffer, MgCl_2 , and Taq enzyme.
Please use aerosol plugged tips to minimize accidental DNA transfer. Be ready to PCR immediately after aliquoting to minimize primer dimers.

• Template (100 ng miniprep DNA):	1 μl
• 5'-primer:	1 μl (1 μM)
• 3'-primer:	1 μl
• Taq mix (Sigma):	25 μl
• MilliQ Water:	<u>22 μl</u>
TOTAL :	50 μl

- 3) Use PCR program 'PCR RRM'. Adjust extension time for the length of your amplified fragment: 1 min. per 1 kb.
 - i. Heated lid: 105°C
 - ii. Preheat lid: on
 - iii. Pause: off
 - iv. Initial denaturation: 94°C 1 min.
 - v. Hot start: off
 - vi. =====Loop 1=====
 - vii. Number of cycles 4
 - viii. Segment 94°C 1 min.
 - ix. Segment 68°C 3 min. 30 sec.
 - x. =====End loop=====
 - xi. =====Loop 2=====
 - xii. Number of cycles 26
 - xiii. Segment 94°C 50 sec.
 - xiv. Segment 62°C 1 min.
 - xv. Segment 70°C 1 min. (adjust)
 - xvi. =====End loop=====
 - xvii. Final extension 72°C 5 min.
 - xviii. Final Hold 4°C