

## **Alcian Blue Hematoxylin/ Orange G Stain for Bone and Cartilage**

**TECHNIQUE:** Formalin fixed, paraffin tissue sections

**REAGENTS:**

**Mayer's Acid Hematoxylin (6 months)**

Hematoxylin (CAS# 517-28-2) -----	2.5 g
Distilled Water -----	350 ml
Aluminum Ammonium Sulfate (CAS# 7784-26-1) -----	25 g
Sodium Iodate (CAS# 7681-55-2) -----	0.25 g
Glycerol (CAS# 56-81-5) -----	150 ml
Glacial Acetic acid (CAS# 64-19-7) -----	10 ml

Dissolve 2.5 g Hematoxylin in 350 ml distilled water with stirring overnight. Next day add 25 g Aluminum Ammonium Sulfate and 0.25 g Sodium Iodate and stir overnight. Next day add glycerol and glacial acetic acid. Filter before use.

**Alcian Blue Hematoxylin (3 months)**

Alcian Blue (CAS# 75881-23-1) -----	1 g
Mayer's Acid Hematoxylin -----	100 ml


**Acid-Alcohol**

70% Ethanol -----	500 ml
Hydrochloric acid, [36.5-38%] -----	5 ml

**0.5% Ammonium Water**

Ammonium Hydroxide, [28.0-30.0%] -----	2.5 ml
Distilled Water -----	500 ml

# CENTER *for* MUSCULOSKELETAL RESEARCH



## **Eosin/Orange G with Phloxine B**

### ***Stock Eosin (2 months)***

Eosin Y (CAS# 17372-87-1) -----	0.6 g
Distilled Water -----	50 ml
100% Ethanol -----	450 ml

Stir well until dissolved. Adjust the pH with glacial acetic acid until the pH is between 4.6 and 5.0

### ***1% Phloxine B (make fresh)***

Phloxine B (CAS# 18472-87-2) -----	0.4 g
Distilled Water -----	40 ml

### ***2% Orange G (make fresh)***

Orange G (CAS# 1936-15-8) -----	0.4 g
Distilled Water -----	20 ml


### ***Working Eosin/Orange G Solution (2 months)***

Stock Eosin -----	500 ml
1% Phloxine B -----	37 ml
2% Orange -----	16 ml



UNIVERSITY of  
**ROCHESTER**  
MEDICAL CENTER

MEDICINE of THE HIGHEST ORDER



## **Alcian Blue Hematoxylin/ Orange G Stain** *(For excellent Growth Plate cartilage morphology)*

**Paraffin sections – 3 microns**

### **PROCEDURE:**


1. Deparaffinize slides and hydrate to distilled water
2. Place in **Acid-alcohol** for **30 seconds** and drain briefly on paper towel (DO NOT RINSE)
3. Place in **Alcian Blue Hematoxylin** for **25 minutes**
4. Wash gently in **distilled water** until excess stain stops leaching from tissue (approximately 3 changes)
5. Differentiate in **Acid-alcohol** for **3 seconds**
6. Rinse gently in **distilled water** with **3 changes**
7. Place in **0.5% Ammonium water** for **15 seconds**
8. Rinse in **distilled water** with **2 changes**
9. Place in **95% EtOH** for **1 minute** (DO NOT RINSE)
10. Place in **Eosin/ Orange G** for **1 min 30 sec**
11. Dehydrate with 3 changes of 95% EtOH and 2 changes of 100% EtOH (1 minute per change)
12. Clear with 3 changes of xylene and coverslip

### **RESULTS:**

Activated osteocytes -----	bright blue pericellular ring
Growth plate -----	pale blue to blue
Cartilage -----	blue/ purple (GAG/ proteoglycan)
Bone -----	orange to red
Erythrocytes -----	bright pink
Soft tissues (muscle, tendon, membranes) --	pink to red
Bone marrow -----	dark blue

\*Colors are more intense in Formic acid decalcified sections compared with EDTA decalcified sections





## Alcian Blue Hematoxylin/ Orange G Stain

*(For staining in both the Growth Plate cartilage and Articular Surface cartilage)*

**Paraffin sections – 3 microns**


### PROCEDURE:

1. Deparaffinize slides and hydrate to distilled water
2. Place in **Acid-alcohol** for **30 seconds** and drain briefly on paper towel (DO NOT RINSE)
3. Place in **Alcian Blue Hematoxylin** for **40 minutes**
4. Wash gently in **distilled water** until excess stain stops leaching from tissue (approximately 3 changes)
5. Differentiate in **Acid-alcohol** for **3 seconds**
6. Rinse gently in **distilled water** with **3 changes**
7. Place in **0.5% Ammonium water** for **15 seconds**
8. Rinse in **distilled water** with **2 changes**
9. Place in **95% EtOH** for **1 minute** (DO NOT RINSE)
10. Place in **Eosin/ Orange G** for **1 min 30 sec**
11. Dehydrate with 3 changes of 95% EtOH and 2 changes of 100% EtOH (1 minute per change)
12. Clear with 3 changes of xylene and coverslip

### RESULTS:

Activated osteocytes -----	bright blue pericellular ring
Growth plate -----	blue
Cartilage -----	blue/ purple (GAG/ proteoglycan)
Bone -----	orange to red
Erythrocytes -----	bright pink
Soft tissues (muscle, tendon, membranes) --	pink to red
Bone marrow -----	dark blue

\*Colors are more intense in Formic acid decalcified sections compared with EDTA decalcified sections



## Alcian Blue Hematoxylin/ Orange G Stain

*(For staining in both the Growth Plate cartilage and Articular Surface cartilage)*

**EMBRYONIC Paraffin sections – 5 or 6 microns**

### PROCEDURE:

1. Deparaffinize slides and hydrate to distilled water
2. Place in **Acid-alcohol** for **20 seconds** and drain briefly on paper towel (DO NOT RINSE)
3. Place in **Alcian Blue Hematoxylin** for **10 minutes**
4. Wash gently in **distilled water** until excess stain stops leaching from tissue (approximately 3 changes)
5. Differentiate in **Acid-alcohol** for **3 seconds**
6. Rinse gently in **distilled water** with **3 changes**
7. Place in **0.5% Ammonium water** for **15 seconds**
8. Rinse in **distilled water** with **2 changes**
9. Place in **95% EtOH** for **1 minute** (DO NOT RINSE)
10. Place in **Eosin/ Orange G** for **1 minute**
11. Dehydrate with 3 changes of 95% EtOH and 2 changes of 100% EtOH (1 minute per change)
12. Clear with 3 changes of xylene and coverslip

### RESULTS:

Activated osteocytes -----	bright blue pericellular ring
Growth plate -----	pale blue to blue
Cartilage -----	blue/ purple (GAG/ proteoglycan)
Bone -----	orange to red
Erythrocytes -----	bright pink
Soft tissues (muscle, tendon, membranes) --	pink to red
Bone marrow -----	dark blue

\*Colors are more intense in Formic acid decalcified sections compared with EDTA decalcified sections