



### Collagen Type 10a1

Quartett Cat # 2031501005.

Immunohistochemistry Protocol For Formalin Fixed Paraffin Embedded Tissue  
Cut sections at 3 microns and bake overnight at 60°C

#### DAY 1

1. Bake slides at **60°C for 30 minutes** prior to starting IHC
2. Deparaffinize slides in xylenes for 5min each, and rehydrate through graded alcohols (100% EtOH twice, 95% EtOH twice and 70% EtOH for 5 min each)
3. Wash twice in deionized water for 5 min each.
4. Enzymatic Antigen Retrieval: Perform antigen retrieval using **Pepsin** (*Sigma P-7000*)
5. (See Pepsin Digest in notes below)
6. Incubate the slides at **37°C** in a water bath **for 10 minutes**.
7. Rinse off pepsin solution in deionized water.
8. Outline each section with a **PAP pen**.
9. Quench endogenous peroxidase in **DAKO Endogenous Blocking Reagent** (*Dako S2003*) **for 30 min**.
10. Rinse thrice with deionized water then once in 1X PBST.
11. Block non-specific binding sites with **the diluted Vectastain Normal Horse Serum** (*Vectastain Elite Universal Kit PK-6200*) – **incubate for 30 min**.
12. **Do not** rinse the serum off after the incubation period.
13. Drain and incubate **overnight at 4°C** with a **1:100 dilution of Collagen X primary antibody** (*Quartett Cat # 2031501005*).
14. Use a **1:200 dilution for embryo sections**. Prepare primary antibody in **2% Normal Horse Serum**. Negative Control slides need to be incubated with **2% Normal Horse Serum** only.



## DAY 2

1. Let slides warm up to room temperature for 30 minutes, then wash 5 times with 1X PBST for 5 min each.
2. Incubate with **Vectastain Biotinylated Horse Anti-Universal secondary antibody** (*Vectastain Elite Universal Kit PK-6200*) for **30 minutes**.
3. Reconstitute the **Vectastain ABC reagent** (*Vectastain Elite Universal Kit PK-6200*) and incubate at **room temperature for 30 min**.
4. Wash 5 times with 1X PBST (5 min each).
5. Incubate the sections with **Vectastain ABC reagent** for **30 minutes**.
6. Wash well with 4 changes of 1X PBST, followed by 2 changes of deionized water
7. Detect color with **Vector Impact DAB** (*Vector SK-4105*) for **a few minutes** (check under the microscope)
8. Stop the color reaction with deionized water.
9. Counterstain the sections with **Hematoxylin** (*Zymed Cat # 93-3943*) for 5 minutes.
10. Rinse with tap water
11. Place slides in 1X PBS for 1-3 minutes.
12. Rinse with deionized water.
13. Dehydrate through 2 changes of 95% ethanol and 2 changes of 100% ethanol
14. Clear in 3 changes of Xylene - mount with Cytoseal.

### ¶ Buffers:-

#### **0.01 N HCl:**

1000 ml of deionized H<sub>2</sub>O + 900  $\mu$ l of concentrated (12N) HCl

#### **Pepsin Digest:**

Dissolve 1 gm of Pepsin powder (*Sigma P-7000*) in 250 ml of 0.01 N HCl solution pre-heated to 37°C. Add the entire solution into a plastic stain bucket. Place the slides in a slide rack and slowly immerse the rack into the Pepsin solution. Cover the stain bucket with its lid and place in 37°C water bath for 10 minutes.

#### **Vectastain Elite Universal IgG Kit (PK-6200):**

**1X PBST** can be used as the diluent buffer for the blocking serum and secondary antibody. Follow the instructions written on the kit data sheet.

¶ **Evidence of positive staining:-**

Collagen type X expression should be localized in the hypertrophic zones of the Growth Plate, as well as the stromal osteoblastic cell populations within the trabecular bone.

Revised on 10/28/2010 by Ashish Thomas, M.S.