



PPAR γ

Cell Signaling # 2435S

Immunohistochemistry Protocol For Formalin Fixed Paraffin Embedded Tissue

Cut sections at 3 microns and bake overnight at 60°C

DAY 1

1. Bake slides at **60°C for minimum of 30 minutes** prior to staining.
2. Deparaffinize tissues in Xylene (3 changes, 5 minutes each), and rehydrate through graded ethanols (2 changes in 100% and 95%, 1 change in 70% for 5 minutes each).
3. Wash twice in deionized water for 5 minutes each.
4. Perform antigen retrieval using a **10mM Sodium Citrate Buffer pH 6.0** – This is done in a water bath at a constant temp. **Set the temp to 65°C** – leave slides in the hot water bath for **1 hour**. Then remove cover and leave outside to cool for about 10- 20 min.
5. Rinse in 3 changes of deionized water.
6. Outline each section with a **PAP pen**.
7. Quench endogenous peroxidase **DAKO Dual Endogenous Enzyme Blocking Reagent for 30 min (Dako S2003)**
8. Rinse with 3 changes of deionized water and once in PBST.
9. Block non-specific binding sites with diluted **Vectastain Normal Goat Serum (Vectastain Elite Rabbit IgG Kit PK-6101) for 30 minutes**.
10. **Do NOT Rinse Slides**. Drain off serum and incubate **overnight at 4°C** with a **1:100 dilution of PPAR γ primary antibody (Cell Signaling # 2435S) for 2-3 month old adult sections**.
11. Prepare the primary antibody in **2% Normal Goat Serum**, this decreases non-specific staining. Negative Control slides need to be incubated with 2% Normal Goat Serum only.



DAY 2

1. Let slides warm up to room temperature for 30 minutes, then wash 5 times with 1X PBST for 5 min each.
2. Incubate with **1:200 Biotinylated Goat anti-rabbit** (*Vectastain Elite Rabbit IgG Kit PK-6101*) secondary antibody for **30 minutes**.
3. Reconstitute the **Vectastain ABC reagent** and incubate at **room temperature for 30 min** (*Vectastain Elite Rabbit IgG Kit PK-6101*)
4. Wash 5 times with 1X PBST for 5 minutes each.
5. Incubate with the **Vectastain ABC reagent** for **30 minutes**.
6. Wash 4 times with 1X PBST, then twice in deionized water.
7. Detect color reaction with **Vector Impact DAB** (*Vector SK-4105*) for a few minutes (check under microscope)
8. Stop the reaction with deionized water.
9. Counterstain the sections with **Hematoxylin** (*Zymed Cat # 93-3943*) for 90 sec.
10. Rinse in tap water for 90 sec.
11. Place slides in 1X PBS for 90 sec.
12. Rinse with deionized water.
13. Dehydrate through 95% ethanol (3 changes) and 2 changes of 100% ethanol.
14. Clear in 3 changes of xylene and mount with cytoaseal.

¶ Buffers:-

10mM Sodium Citrate, 0.05% Tween-20 Buffer; pH 6.0

1. Tri-Sodium Citrate Dihydrate 2.94g
2. dH₂O 1000 ml

Mix well to dissolve. Adjust pH to 6.0 with 1N HCl

Add 0.5ml of Tween 20 and mix. Store the buffer at 4°C for longer storage.

¶ **Vectastain Elite Rabbit IgG Kit (PK-6101):**

1XPBST can be used as the diluent buffer for the blocking serum and secondary antibody. Follow the instructions written on the kit data sheet.