



## Anti -PTH/ PTHrP Receptor

(Upstate Cell Signaling Solutions Catalog # 05-517 - Mouse monoclonal antibody)

### **Immunohistochemistry Protocol For Formalin Fixed Paraffin Embedded Tissue Cut sections at 3 microns and bake overnight at 60°C**

#### **DAY 1**

1. Bake slides at **60°C for minimum of 30 minutes** prior to staining.
2. Deparaffinize tissues in xylene (3 changes for 5 minutes each), and rehydrate through graded ethanols (2 changes of 100% and 95%, 1 change 70% for 5 minutes each).
3. Wash twice in deionized water for 5 minutes each.
4. Perform antigen retrieval using a **10mM Sodium Citrate Buffer pH 6.0** – This is done in a water bath at a constant temp. **Set the temp to 65°C** – leave slides in the hot water bath for **1 hour**. Then remove cover and leave outside to cool for about 10- 15 min.
5. Rinse in 3 changes of deionized water for 3-5 minutes each.
6. Outline each section with a **PAP pen**.
7. Quench endogenous peroxidase **DAKO Dual Endogenous Enzyme Blocking Reagent for 30 min (Dako S2003)**
8. Rinse briefly with 2 changes of PBST.
9. Block endogenous **Avidin** with the **Vector Avidin-Biotin Blocking Kit (Vector SP-2001) for 15 minutes**.
10. Rinse briefly with 2 changes of PBST.
11. Block endogenous **Biotin** with the **Vector Avidin-Biotin Blocking Kit (Vector SP-2001) for 15 minutes**.
12. Rinse briefly with 2 changes of PBST.
13. Block non-specific binding sites with the working solution of **M.O.M Mouse Ig Blocking Reagent (Vector M.O.M kit # BMK-2202)** prepared as described in the data sheet, for **1 hour**.
14. Wash sections twice in PBST for 2 minutes each.
15. Incubate sections in the working solution of **M.O.M diluent (Vector M.O.M kit # BMK-2202) for 5 minutes**.
16. **Do not Rinse Slides!!** Drain off M.O.M diluent and incubate **overnight at 4°C** with a **1:100 dilution of Anti-PTH/PTHrp primary antibody (Upstate Cell Signaling Solutions Cat # 05-517)**.
17. . Negative Control slides need to be incubated with the M.O.M diluent only.



## DAY 2

1. Let slides warm up to room temperature for 15-20 minutes, then wash twice with PBST for 5 min each.
2. Incubate with the working solution of **M.O.M Anti-Mouse Secondary Antibody** (*Vector M.O.M kit # BMK-2202*) for **10 minutes**.
3. Wash twice with PBST for 2 minutes each.
4. Reconstitute the **Vectastain ABC reagent** and incubate at **room temperature for 30 min** (*Vectastain Elite Mouse IgG Kit PK-6102*)
5. Incubate with the **Vectastain ABC reagent** for **30 minutes**.
6. Wash twice with PBST, then twice in deionized water for 5 minutes each.
7. Detect color reaction with **Vector Impact DAB** (*Vector SK-4105*) for a few minutes (check under microscope)
8. Stop the reaction with deionized water.
9. Counterstain the sections with **Hematoxylin** (*Zymed Cat # 93-3943*) for 10 minutes.
10. Wash in tap water for 5 minutes.
11. Place slides in 1X PBS for 1-3 minutes.
12. Rinse with deionized water.
13. Dehydrate through 95% ethanol (3 changes) and 2 changes of 100% ethanol.
14. Clear in 3 changes of xylene and mount with cyto seal.

### Buffers:-

#### **10mM Sodium Citrate Buffer; pH 6.0**

1. Tri-Sodium Citrate Dihydrate      2.94g
2. dH<sub>2</sub>O                                      1000 ml

Mix well to dissolve. Adjust pH to 6.0 with 1N HCl

Add 0.5ml of Tween 20 and mix. Store the buffer at 4°C for longer storage.

### Evidence of Positive Staining.

Expression should be localized to the pre-hypertrophic zones of the growth plate, early hypertrophic chondrocytes, stromal osteoblast cell populations in the trabecular bone and in the deeper zones of the articular cartilage. It should **not** be expressed in the columnar cells of the growth plate.

**(Revised on 10/17/2011 by Ashish Thomas)**