

## TCID 50 protocol

1. Put A6 cells in 96well plates at 5000cells/well in ASF +10% FCS (100ul total) for about 2 days.
2. In a separate V bottom 96 well plate, do dilutions of virus:
  - A. 1:25- Add 8ul of virus in 192ul of ASF-A6 or ASF
  - B. 1:50- Add 60ul of 1:25 in 60ul ASF-A6 or ASF
  - C. 1:500- Add 12ul of 1:50 in 108ul ASF-A6 or ASF
  - D. 1:5000- Add 12ul of 1:500 in 108ul media and so on until dilution H.  $5 \times 10^7$

\*Note: Make sure you resuspend the wells very well before going to the next dilution!!
3. Transfer 100ul of each dilution to the plate with A6 to make a total of 200ul per well.

\*Remember to leave columns 1 and 2 without virus, these are the control.
4. Leave in incubator at 37C with CO2 for approx. 5 days.
5. Check daily.

## TCID 50 Calculation

### Example

<u>Dilution</u>		<u>Control</u>												<u>% Death</u>
		1	2	3	4	5	6	7	8	9	10	11	12	
$10^{-1}$	A	(-)	(-)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	100%
$10^{-2}$	B	(-)	(-)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	100%
$10^{-3}$	C	(-)	(-)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	100%
$10^{-4}$	D	(-)	(-)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	100%
$10^{-5}$	E	(-)	(-)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	100%
$10^{-6}$	F	(-)	(-)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	(+)	(-)	(-)	80%
$10^{-7}$	G	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	0%
$10^{-8}$	H	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	(-)	0%

**1. Calculate Proportionate Distance (PD) between the two dilutions in between 50% death: (% next above 50%)- 50% / (% next above 50%) – (% next below 50%)**

Example above:

$$PD = 80\% - 50\% / 80\% - 0\% = 30/80 = .375$$

**2. Calculate 50 % end point. Log lower dilution= dilution in which position is next above 50%**

Example above:

$$\text{Log lower} = 10^{-6} \text{ or } -6$$

**3. Add PD and Log lower dilution**

$$\text{Example above: } -6 + .375 = -6.375$$

$$\text{Log TCID}_{50} = 10^{-6.375} \text{ or } 1 / 2.37 \times 10^6$$

**4. Calculate TCID 50/ml. Divide by the ml of viral inoculum added to row A**

$$\text{Example above: according to our protocol} = .008 \text{ ml}$$

$$\text{TCID}_{50}/\text{ml} = 2.37 \times 10^6 / .008 = 2.9 \times 10^8$$

**5. Calculate PFU/ml. Divide by constant.**

$$\text{Example above: } 2.37 \times 10^6 \times .69 = 2.0 \times 10^8 \text{ PFU/ml}$$