

Functional Genomics Center Bioanalyzer Only Sample Drop-off Form

Please note this form is ONLY used for samples the FGC will run on the Agilent Bioanalyzer

Investigator (PI):	Date:
Contact Name:	Phone#:
Account#	Reference# (if needed for billing)
Email:	
Billing Contact Information:	

(circle one) Prokaryotic / Eukaryotic

Please submit at least 3 uL of sample; 5ul if you don't know the concentration (ng/ul)

	Sample Name	Sample storage buffer (ddH2O, buffer)	[sample] ng/ul	Volume ul	RNA	Small RNA	DNA	If DNA, size range?
1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								
13								
14								
15								
16								
17								
18								
19								
20								

***Definitions:**

RNA -Based on the concentration, we will choose to run samples with either the RNA 6000 Nano Kit (25-500ng/ul) or the RNA 6000 Pico Kit (0.05-25ng/ul)

Small -Small RNA Kit; detailed view of 6-150 nt sizes; concentrations of 1-100ng/ul for total RNA and 1-20ng/ul for enriched small RNA

DNA -Based on the concentration and size range, we will choose to run samples with either the High Sensitivity DNA kit (50-7000 bp; 5-500pg/ul), DNA1000 (25-1000bp; 0.1-50ng/ul), DNA12000 (100-12000bp; 0.5-50ng/ul) Different fees may apply;

Questions? call 276-9981 or email corelab@urmc.rochester.edu

For Core Use Only	Comments:
Date Run: Initials:	
Email Results Date: Initials:	
Billing Date: Initials:	