

Image Exploration

We recommend the use of ImageJ, a versatile freeware program for medical imaging analysis. It is able to run both on Macs, Windows and Linux operating systems and is able to import the Fluoview Native File formats(.oib and .oif) with the installation of a plugin.

To Set-Up ImageJ:

- 1- Go to <http://rsbweb.nih.gov/ij/>
- 2- Click on the download link
- 3- Download the appropriate version for your machine
- 4- Install Software into default directory
- 5- Go to <http://www.loci.wisc.edu/ome/formats-imagej.html>
- 6- Click the download loci-tools.jar hyperlink
- 7- Place the downloaded .jar file into the plugins folder of ImageJ

Documentation:

Instructions for the use of ImageJ can be found online at <http://rsbweb.nih.gov/ij/docs/index.html> including in depth tutorials