

New Shelf Life for CLIDOX™ (Chlorine Dioxide) Solution

Our vendor of Clidox™ informed us that recent studies indicate their high level disinfectant concentration of chlorine dioxide (1:5:1) loses no efficacy over two weeks. Before January 1, 2009, the vivarium will implement a new schedule for making up fresh Clidox™ in rodent housing and procedure rooms documenting a two week shelf life of the reconstituted solution INSTEAD OF 24 HOURS. Please don't hesitate to contact me if I may provide more information.

Jeff Wyatt DVM, MPH, DACLAM
Professor & Chair
Department of Comparative Medicine
University of Rochester Box 674
Rochester, New York 14642

T: 585-275-2651

F: 585-273-1085

E: jeff_wyatt@urmc.rochester.edu