

Thromboelastograph Testing (TEG)

To: Clients of UR Medicine Labs

From: Majed Refaai, M.D.
Associate Professor of Pathology & Laboratory Medicine
Special Coagulation Laboratory Director
University of Rochester Medical Center
Phone: (585) 276-3927 email: Majed_Refaai@urmc.rochester.edu

Robert J. Miller
Chief Supervisor, UR Medicine, Strong Memorial Hospital, Hematology &
Chemistry Laboratory
Phone: (585) 275-3244 email: Robert_Miller@urmc.rochester.edu

Date: July 20, 2015

The TEG Thromboelastograph test has been added to the test menu of the UR Medicine Labs Hematology and Chemistry Laboratory. The TEG analyzer uses whole blood to measure all phases of hemostasis. There are two orderable tests: TEG and HTEG (Heparinized TEG). HTEG should be used for patients on heparin.

Sample requirement: 1 full NaCitate (Blue Top) tube.

Test availability: Available at all times, 24/7.

Sample Stability Must be performed within 2 hrs. of collection.

Go to www.testmenu.com/Rochester for additional collection guidelines.