

## Celiac Disease Screen: Change to Assay

Date: November 7, 2017

---

Effective Date: December 4, 2017

The new reference range for the Celiac IgA assay is:

<15 U/mL – NEGATIVE  
≥15 U/mL – POSITIVE

Note: Total IgA will no longer be included in the Celiac screen. If the assay suggests that the patient is IgA-deficient, tissue Transglutaminase (tTG) IgG and deaminated gliadin peptide (DGP) IgG will be reflex ordered as an alternative screen for Celiac disease.

From: Dwight Hardy, PhD  
Director, Microbiology Laboratory, SMH  
Phone: (585) 275-1408  
Email: [Dwight\\_Hardy@urmc.rochester.edu](mailto:Dwight_Hardy@urmc.rochester.edu)