

ELECTRODE PRODUCT GUIDE








Product Name	Differentiation Features	Location Used	Defibrillator
CPR-D-padz 	<ul style="list-style-type: none"> Enables Real CPR Help® and See-Thru CPR® One-piece design simplifies proper positioning on the chest One size fits all, with ability to expand for larger chest 	<ul style="list-style-type: none"> EMS & Fire Public Access 	<ul style="list-style-type: none"> M Series® E Series® R Series® AED Plus® AED Pro® X Series®
CPR Stat-padz 	<ul style="list-style-type: none"> Enables Real CPR Help and See-Thru CPR Multi-function, two-piece design ensures all therapies can be delivered Accommodates AP and AA pad placement 	<ul style="list-style-type: none"> Hospital Med/Surg floors ER & ICU EMS & Fire 	<ul style="list-style-type: none"> M Series E Series R Series AED Plus AED Pro X Series
OneStep Basic 	<ul style="list-style-type: none"> Multi-function electrode Automated self-test with R Series 	<ul style="list-style-type: none"> Hospital 	<ul style="list-style-type: none"> M Series R Series X Series
OneStep Pacing 	<ul style="list-style-type: none"> OneStep Basic features plus: <ul style="list-style-type: none"> OneStep Pacing 	<ul style="list-style-type: none"> Hospital 	<ul style="list-style-type: none"> R Series
OneStep CPR (AP) 	<ul style="list-style-type: none"> OneStep Basic features plus: <ul style="list-style-type: none"> Real CPR Help See-Thru CPR Anterior-Posterior placement 	<ul style="list-style-type: none"> Hospital 	<ul style="list-style-type: none"> M Series R Series X Series
OneStep CPR (AA) 	<ul style="list-style-type: none"> OneStep Basic features plus: <ul style="list-style-type: none"> Real CPR Help See-Thru CPR Anterior-Anterior placement 	<ul style="list-style-type: none"> Hospital 	<ul style="list-style-type: none"> M Series R Series X Series
OneStep Complete 	<ul style="list-style-type: none"> All OneStep Basic, Pacing, and CPR capabilities 	<ul style="list-style-type: none"> Hospital 	<ul style="list-style-type: none"> R Series
OneStep Pediatric 	<ul style="list-style-type: none"> Designed for patients 0–8 years of age and less than 55 pounds (25 kg) Automatically adjusts joules setting to pediatric levels 	<ul style="list-style-type: none"> Hospital 	<ul style="list-style-type: none"> M Series R Series X Series
OneStep Pediatric CPR 	<ul style="list-style-type: none"> Designed for patients 0–8 years of age and less than 55 pounds (25 kg) Automatically adjusts joules setting to pediatric levels Real CPR Help See-Thru CPR 	<ul style="list-style-type: none"> Hospital 	<ul style="list-style-type: none"> M Series R Series X Series
Pedi-padz® II 	<ul style="list-style-type: none"> Connects directly to AEDs & ALS devices Pediatric intelligence function with ZOLL AEDs 	<ul style="list-style-type: none"> EMS & Fire Public Access 	<ul style="list-style-type: none"> M Series E Series R Series AED Plus AED Pro X Series
Pedi-padz Multi-function Liquid or Solid Gel 	<ul style="list-style-type: none"> Multi-functional electrode for pediatric patients Available in liquid or solid gel 	<ul style="list-style-type: none"> Hospital EMS & Fire 	<ul style="list-style-type: none"> M Series E Series R Series X Series
Pedi-padz Radiolucent Solid Gel 	<ul style="list-style-type: none"> For use in cath and RP labs for procedures Virtually invisible to x-ray 	<ul style="list-style-type: none"> Hospital <ul style="list-style-type: none"> Cardiac cath lab Electrophysiology labs 	<ul style="list-style-type: none"> M Series E Series R Series X Series
Pro-padz® Multi-function Liquid or Solid Gel 	<ul style="list-style-type: none"> Available in liquid or solid gel Available with or without tray Used in cath and EP labs for procedures 	<ul style="list-style-type: none"> Hospital <ul style="list-style-type: none"> Cardiac catheterization lab Electrophysiology labs 	<ul style="list-style-type: none"> M Series E Series R Series X Series
Pro-padz Radiolucent Liquid or Solid Gel 	<ul style="list-style-type: none"> For use in cath and EP labs for procedures - Virtually invisible to x-ray Pre-Connect (solid only) Speed pack 	<ul style="list-style-type: none"> Hospital <ul style="list-style-type: none"> Cardiac catheterization lab Electrophysiology labs 	<ul style="list-style-type: none"> M Series E Series R Series X Series
Pro-padz Sterile 	<ul style="list-style-type: none"> For use in operating rooms Helps maintain sterile barrier Packaged in two pouches Available in 54" or 10-foot lead wires 	<ul style="list-style-type: none"> Hospital <ul style="list-style-type: none"> Operating rooms 	<ul style="list-style-type: none"> M Series E Series R Series X Series
Stat-padz® 	<ul style="list-style-type: none"> Multi-function Pre-Connect Speed pack Different pad shapes and colors 	<ul style="list-style-type: none"> EMS & Fire Hospital 	<ul style="list-style-type: none"> M Series E Series R Series X Series
Stat-padz II 	<ul style="list-style-type: none"> Connects directly to AEDs and ALS devices 	<ul style="list-style-type: none"> EMS & Fire Hospital 	<ul style="list-style-type: none"> M Series E Series R Series AED Plus AED Pro X Series

See reverse side for proper cable/connector requirements.

ELECTRODE GUIDE



	CPR-D-padz®	CPR Stat-padz®	OneStep™ Complete	OneStep™ CPR AP	OneStep™ CPR AA	OneStep™ Pediatric CPR
AED Plus®	 <ul style="list-style-type: none"> Connects directly for Real CPR Help® and defibrillation 	<ul style="list-style-type: none"> Connects directly for Real CPR Help and defibrillation 	<ul style="list-style-type: none"> Not compatible 	<ul style="list-style-type: none"> Not compatible 	<ul style="list-style-type: none"> Not compatible 	<ul style="list-style-type: none"> Not compatible
AED Pro®	 <ul style="list-style-type: none"> Connects directly for Real CPR Help, See-Thru CPR®, and defibrillation 	<ul style="list-style-type: none"> Connects directly for Real CPR Help, See-Thru CPR, and defibrillation 	<ul style="list-style-type: none"> Not compatible 	<ul style="list-style-type: none"> Not compatible 	<ul style="list-style-type: none"> Not compatible 	<ul style="list-style-type: none"> Not compatible
X Series®	 <ul style="list-style-type: none"> Requires CPR Connector (PN# 8000-0370) for CPR feedback Connects directly for defibrillation 	<ul style="list-style-type: none"> Requires CPR Connector (PN# 8000-0370) for CPR feedback Connects directly for defibrillation 	<ul style="list-style-type: none"> Connects directly for all capabilities No OneStep Pacing Capabilities 	<ul style="list-style-type: none"> Connects directly for all capabilities 	<ul style="list-style-type: none"> Connects directly for all capabilities 	<ul style="list-style-type: none"> Connects directly for all capabilities
E Series®	 <ul style="list-style-type: none"> Requires CPR Connector (PN# 8000-0370) for CPR Feedback Connects directly for defibrillation 	<ul style="list-style-type: none"> Requires CPR Connector (PN# 8000-0370) for CPR feedback Connects directly for defibrillation 	<ul style="list-style-type: none"> Not compatible 	<ul style="list-style-type: none"> Not compatible 	<ul style="list-style-type: none"> Not compatible 	<ul style="list-style-type: none"> Not compatible
R Series®	 <ul style="list-style-type: none"> Requires R Series® CPR Connector (PN# 8009-0020) for CPR Feedback Connects directly for defibrillation 	<ul style="list-style-type: none"> Requires R Series CPR Connector (PN# 8009-0020) for CPR Feedback Connects directly for defibrillation 	<ul style="list-style-type: none"> Connects directly for all capabilities 	<ul style="list-style-type: none"> Connects directly for all capabilities 	<ul style="list-style-type: none"> Connects directly for all capabilities 	<ul style="list-style-type: none"> Connects directly for all capabilities
Shelf Life	5 yrs	2 yrs	2 yrs	2 yrs	2 yrs	2 yrs
Quantity and Part Numbers	Individual: 8900-0800-01	Individual: 8900-0402 8 per case: 8900-0400	Individual: 8900-0224-01 8 per case: 8900-0214-01	Individual: 8900-0223-01 8 per case: 8900-0213-01	Individual: 8900-0225-01 8 per case: 8900-0217-01	Individual: 8900-000219-01 8 per case: 8900-000220-01

(AP: anterior-posterior; AA: anterior-anterior)
See reverse side for individual product details.