

September, 2016

Vitamin D Testing

As of August 26, 2016 testing of Vitamin D by LCMS/MS returned to Strong Memorial Hospital Clinical Toxicology/Mass Spectrometry Laboratory.

The reference interval for the SMH Clinical Toxicology/Mass Spectrometry Vitamin D test is as follows:

SMH Reference Interval: 25-Hydroxyvitamin D₂ and D₃ by Tandem Mass Spectrometry, Serum

Deficiency	Less than 20 ng/mL
Low risk of deficiency	20-29 ng/mL
Adequate	30-60 ng/mL
Potentially harmful	Greater than 60 ng/mL

Communication

- We are updating e-mail addresses in our database. If you have an additional e-mail you would like us to have, please send it to : jessica_wahl@urmc.rochester.edu

Patient Service Center (PSC) Updates

- **RELOCATION: PSC at 2030 Monroe Ave. has moved to 1520 Monroe Ave.**
This location offers improved parking and first floor entry.
The hours remain the same M-F 8am-5pm, (closed 12pm-1pm)
- **NEW PSC in LeRoy:** 127 W. Main St., Suite 5 LeRoy NY, 14482 (across the street from Tops)
Hours: M-F 7:30am-4:30pm (closed 12pm-1pm)
- **NEW PSC in Downtown Rochester:** 454 East Broad St., Suite 200, Rochester, NY 14607
Hours: M-F 7am-4pm (closed 12pm-1pm)

How to contact us: email: labservicesoutreach@urmc.rochester.edu

Client Services Call Center: (585) 758-0510, option 3

Go to: www.urmedicinelabs.org

Test Menu: www.testmenu.com/rochester